


How to Export Scans as .STL Files from a Lythos Scanner

Log into your OrmcoDigital Account

Access your scans with a click of a button by visiting www.ormcodigital.com

If you are a new user or have never accessed ormcodigital.com before, please contact Ormco Customer Support: **800.854.1741 (option 2; option 1)** to receive access to your account or setup a new one.



Patient List Window

The first window you will see is the “**Patient List**”, which shows users all cases that have been submitted for an Ormco Custom product. The tabs on this window include valuable information:

- **ALERTS:** Shows all cases that have been placed on hold
- **PENDING APPROVAL:** shows all cases that require Doctor approval prior to proceeding
- **ALL:** shows all cases that are currently being manufactured or designed by Ormco
- **HISTORY:** shows all cases that have been completed

Welcome, Dr. Units




Image	Patient Name	Case ID	Case Date	Bonding Date	Current Status	Approver
Image Not Available	Michael Panthin	38161	9/6/13	10/2/13	Awaiting Approval	

You can track the progress of your case by using the search option in the dropdown menu and entering the keyword or number you'd like to search for.




Lythos Scans Window

To access new Lythos scans click on the “Lythos Scans” tab located in the middle of the screen:



This window includes all scans that have been uploaded from your Lythos scanner:



If you notice that a case isn't showing up in this window, please make sure that your case uploaded successfully on your Lythos scanner. If the case failed, resend the impression. Once it has uploaded, you should see the case move from the “In-process case transfers” window to the “Completed case transfers” window with a “Transferred” result.

If your case has uploaded but isn't showing up on **ormcodigital.com**, then your Lythos device may not be matched properly. Perform the following steps:

1. Log in as a User with Administrator rights.
2. On your Lythos scanner please go to the System Menu by pressing the gear in the top left-hand corner, select **"Administration"** and then select **"Users"**.
3. Click on the Doctor in the list and select **"Edit User"**.
4. Enter the correct Ormco Account Number in the Ormco ID file
5. Save the User Profile
6. Once back in the Add/Edit Users window, press **"Save"** again
7. Now try a simple scan, upload the case, and ensure that it is accessible on your **ormcodigital.com** account

If you are still experiencing problems, please call **Lythos Tech Support at 800-854-1741 (option 2; option 2)**.

Lythos is an Open Platform, which means that you have complete ownership of all completed scans. There are 3 buttons users should use to access scans:

1. **VIEW SCAN:** displays a simple thumbnail image of your case
2. **DOWNLOAD:** allows customers to download the raw scan files in .stl format to your hard drive

Welcome, Dr. Units


Patient List		Lythos Scans					
		Delete Scans		From: <input type="text"/>	To: <input type="text"/>	Search <input type="text"/>	
New Scans 55		Patient Name	Scan Number	Scan Date	View Scan	Orders	Hi-Res Scans
All Scans 61	<input type="checkbox"/>	H.G	DIS137015N-CASE_00014	2/28/14 5:52 PM		No Orders	Download Place Order
<small>Please note: scans with a red flag next to the scan number have been marked as incomplete.</small>							
	<input type="checkbox"/>	J.H	DIS132-PP18-CASE_00112	10/11/13 6:44 PM		No Orders	Download Place Order
	<input type="checkbox"/>	Jeanne Huribout	DIS132-PP18-CASE_00108	10/11/13 5:46 PM		No Orders	Download Place Order
	<input type="checkbox"/>	Merrill Jeanne	DIS133-PP26-CASE_00107	12/17/13 12:06 AM		No Orders	Download Place Order
	<input type="checkbox"/>	H.J	DIS132-PP18-CASE_00113	10/11/13 6:58 PM		No Orders	Download Place Order

Downloading your case

When you click on the **"Download"** button you will receive the following prompt:



Clicking “Download” will initiate the download of a zip file that includes 2 files named “Lower-Arch.stl” and “Upper_Arch.stl”. These files are typically about 15 megabytes per files but are downloaded as a zipped file that is close to 15 megabytes total. Once the files has completed downloading, you can access it from the Downloads folder on your computer:



To access the scans there are two options:

1. Open the zipped folder, copy the two files within the folder, and paste them into another folder
2. Unzip the folder by right-clicking on the folder, selecting “Extract All”, and then selecting a destination for the files

FINAL STEP:

Log on to specialtyappliances.com and upload your files to the Specialty Appliances Digital Portal. If you need assistance, contact our digital customer service at 800-522-4636 or 678-513-4408.