

You and your patient have worked hard to produce a beautiful smile. Despite the best intentions of the most diligent patient, anterior teeth can relapse with removable retainers. Being able to tell your patients they will experience little if any anterior shifting is the primary benefit of using a Fixed Lingual Retainer (FLR) from Specialty Appliances.



Our most popular FLR design utilizes a solid stainless steel wire contoured to the lingual surface of each anterior tooth. The wire is virtually impossible to break and the adhesive pads are kept thin for patient comfort.



Fixed Lingual retainers can also be used on the upper arch. Our typical upper design has a composite pads on the centrals and laterals with the wire carefully bent to avoid anterior interferences.



An often prescribed lower FLR design features a solid round wire contoured to the anterior teeth and bonded to only the cuspids. The patient is able to easily use a floss threader to keep the Retainer clean.



Ideally, FLR's from Specialty are delivered at the same appointment the braces are removed, this reduces the chances of undesirable tooth movement. We can fabricate your FLR's from PVS impressions, stone models or digital scans. This allows immediate delivery when the braces are removed.

"I've been using Fixed Lingual Retainers from Specialty Appliances for over 20 years with a 98% success rate. Specialty has developed a flawless product making the clinical procedure one of the most efficient in my office."

-Dr. Lynn Davis, TX

"Indirect bonding lower cuspid to cuspid retainers has changed the way I practice. It takes little of my time to place it, and the fit is always perfect. I look forward to seeing my retention patients because I know their lower incisors will remain straight."

-Dr. David Sain Murfreesboro, TX

Please visit our website today at www.specialtyappliances.com for additional information on our FLR service.