

Invisible Retainers from Specialty Appliances meet your patients' night time retention needs, at a great value for your practice. We use the industry leading Essix brand material for its exceptional clarity and durability. Take advantage of our new Invisible Retainer Program today!



## Benefits with Specialty Invisible Retainers

- We use the “Industry Standard” Essix material with exceptional clarity and durability
- “Value Pricing” when multiple sets of Invisibles are ordered
- Minor anterior corrections can be added upon request and built into the Invisible Retainers
- For your convenience we accept all intraoral scans for Invisible Retainers

## Here's how the Invisible Retainer program works

*Send us a set of stone models, PVS impressions, or digital impressions from any intraoral scanner. A 3D printed replica of your patients model will be provided at no cost for orders of two or more retainers per arch. We will then ship the retainers to your office. You then have the option of delivering them all at once or reserving extra retainers for future delivery at a “Retainer Check” appointment.*

*A popular retention combination utilizes our Invisible Retainers fabricated to fit over a Fixed Lingual Retainer (FLR). The Invisible Retainers are processed directly over the FLR's to guarantee an accurate fit and maximum patient comfort. The FLR's come with a transfer tray for easy delivery. Utilize this service for either or both arches.*

Please visit our website today at [www.specialtyappliances.com](http://www.specialtyappliances.com) for further information regarding our Invisible Retainer service.