

Splints

Models/Scans: Upper and lower models poured in Class I or Class II stone (models). Complete anterior and posterior anatomy of **both arches** so the appliance can be equilibrated in the lab.

Wax Bites: Wax bites should reflect the exact vertical dimension of interocclusal acrylic in the final appliance. Digital bites can be sent as well.

Notes: Indicate articulator settings when mounted models are sent to the lab. Be careful to watch for severe undercuts. It may be necessary to use Variflex® material which is forgiving in patients with severe undercuts.

