


The M4 Herbst Telescope

Specialty Appliances and the new M4 Herbst Telescope has quickly built a reputation of durability and patient comfort. With 40 degrees of lateral movement and a maximum incisal opening of 65mm, the M4 has almost eliminated restricted movement for the Herbst patient. This has not only improved patient comfort, but has drastically reduced breakage. The M4 telescope is offered in 16mm and 19mm, improving comfort and allowing doctors to bond upper 5-5 and lower 4-4. All anchorage systems work with M4 including bands, crowns and acrylic appliances.



The Difference: M4 allows the patient a maximum incisal opening of 65mm compared to the competitor's 38mm.

To further improve Herbst treatment with the M4 telescope, we recommend the Specialty Appliances AppleCore Screws. We have modified the geometry of the Herbst screw body to incorporate a reduced radius towards the center which allows greater lateral movement of the lower jaw. Simply request AppleCore Screws for any Herbst design and we will provide them at NO ADDITIONAL CHARGE!



Applecore Screws increase lateral movement for normal lateral function