

## Removable Retainers

**Models/Scans:** Upper and lower models poured in Class I or Class II stone (models). Complete anterior and posterior anatomy of the arch on which the appliance will be fabricated. Palatal anatomy is required for the fabrication of upper Hawley retainers.

**\*Opposing Models:** It is strongly suggested that the opposing arch be sent as well when requesting appliances in which bite planes, habit cribs or pontics are added.

**Wax Bites:** Necessary only if centric occlusion is not obvious when articulating models.

**Notes:** Reference the Specialty Appliances color guide when ordering retainers and be sure to include the desired color code on your prescription sheet (Hawley retainers).

