

SEND: Rx FORMS MAILING LABELS SHIPPING SUPPLIES

DOCTOR _____ ACCT# _____

ADDRESS _____

CITY _____ STATE _____ ZIP _____

PHONE _____ EMAIL _____

PATIENT NAME _____

DATE SHIPPED _____ DATE NEEDED* _____

**Date needed should be at least 1 day before appointment date.*

APPROVAL TO CHARGE EXPEDITED MANUFACTURING & SHIPPING

DIGITAL SCAN TAKEN WITH:

- iTero® Carestream 3M Medit TRIOS® Sirona
 Other _____

CLEAR IMAGE® ALIGNERS

Opposing Model Required for Warranty

- | | | |
|---|-----------------------|-----------------------|
| | U | L |
| <input type="radio"/> Select 1 Aligner | <input type="radio"/> | <input type="radio"/> |
| <input type="radio"/> Select 3 Aligners | <input type="radio"/> | <input type="radio"/> |
| <input type="radio"/> Select 5 Aligners | <input type="radio"/> | <input type="radio"/> |
| <input type="radio"/> Plus 10 Aligners | <input type="radio"/> | <input type="radio"/> |
| <input type="radio"/> Refinement Aligners | <input type="radio"/> | <input type="radio"/> |
- (this case has been submitted previously)*

Send Tx Review*

- Yes* Email _____
*By selecting this you will be responsible for a setup fee if patient does not proceed with treatment.
- No
- Add Digital Pontics
 Location _____
 Shade (Paint) _____

CIRCLE TEETH TO RESET

R	7	6	5	4	3	2	1	1	2	3	4	5	6	7	L
	7	6	5	4	3	2	1	1	2	3	4	5	6	7	

INVISIBLE RETAINERS

- | | | | | | |
|---------------------------|-----------------------|-----------------------|-------------|-----------------------|-----------------------|
| | U | L | | U | L |
| Single Invisible Retainer | <input type="radio"/> | <input type="radio"/> | Guardian® 3 | <input type="radio"/> | <input type="radio"/> |
| Guardian® 2 | <input type="radio"/> | <input type="radio"/> | Guardian® 4 | <input type="radio"/> | <input type="radio"/> |

INVISIBLE RETAINER OPTIONS

- Add Fixed Lingual Retainers *Indicate in section* →
- Add Pontics Location _____ Shade _____
- Add Mouthguard Upper Lower
 Clear Black White Blue Red Green

SPECIAL INSTRUCTIONS

DOCTOR SIGNATURE _____


License # _____ Expiration _____

OFFICE USE: 1 2 3 4 + PD: SA DR

MODELS: U L BOTH BANDS CROWNS BROKEN

IMPRESSIONS: U L BOTH

DISINFECT: _____ QA IN: _____ FINAL INSP: _____



IPR INSTRUCTIONS

- Reduce Teeth as Needed in Laboratory
 No IPR Required

SPACE & OVERJET INSTRUCTIONS

- Okay to leave space on upper due to overjet restrictions. Leave space: _____
- Opposing arch unavailable. Overjet measurement is: _____

FIXED LINGUAL RETAINERS (FLR)

- | | | |
|-------------------------------------|-----------------------|-----------------------|
| Placement of Retainer | U | L |
| Central - Central | <input type="radio"/> | <input type="radio"/> |
| Lateral - Lateral | <input type="radio"/> | <input type="radio"/> |
| Cuspid - Cuspid | <input type="radio"/> | <input type="radio"/> |
| Placement of Pads | | |
| Composite Pads on Each Tooth | <input type="radio"/> | <input type="radio"/> |
| Composite Pads on Distal Most Teeth | <input type="radio"/> | <input type="radio"/> |
| Mesh Pads on Each Tooth | <input type="radio"/> | <input type="radio"/> |
| Mesh Pads on Distal Most Teeth | <input type="radio"/> | <input type="radio"/> |
| Type of Wire | | |
| Round .028 | <input type="radio"/> | <input type="radio"/> |
| .016 x .022 Braided | <input type="radio"/> | <input type="radio"/> |
| .016 x .022 Solid Stainless Steel | <input type="radio"/> | <input type="radio"/> |

(800) 522-4636
 FAX (470) 239-7217
 www.specialtyappliances.com

CLEAR IMAGE® ALIGNERS

Select 1, 3, or 5 Aligner(s)

Select Aligners are used for the simplest aligner cases that require minor tooth movements. This product is offered in quantities of 1, 3, and 5 depending on how much tooth movement is required. Select aligners can be ordered for single or dual arch.

Plus 10 Aligners

Plus Aligners are used when moderate tooth corrections are required. Ten aligners are included with this option and can be ordered for single or dual arch.

Send Tx Review

A Tx review will allow for the setup to be sent to the doctor for approval prior to aligner fabrication. If a Tx Review is requested and the patient does not proceed with treatment, a setup fee will be incurred. If a Tx Review is not requested, or neither option is selected, fabrication of the aligners will begin immediately.

Circle Teeth to be Reset

The diagram provided allows the clinician to select individual or groups of teeth for corrections by marking the diagram accordingly.

Opposing Models

Specialty Appliances requires an opposing scan or model be submitted when making Clear Image® Aligners. This ensures no occlusion interference is created by tooth movements and guides us on necessary tooth movements for best outcomes. If an opposing scan or model is not submitted, the warranty for success of Clear Image® Aligners is not offered and any refinement aligners will be subject to full fees.

IPR Instructions — Use Diagram

Most Clear Image® Aligners require IPR to help create space for the teeth to move. In this section, please specify whether Specialty Appliances is to determine IPR amounts and locations or that no IPR is required.

Scanning Stone Models

Clear Image® Aligner cases that are sent in as impressions or stone models will be subject to automatic impression pouring and scanning fees. If the patient chooses not to proceed with treatment, these fees will still apply.

GUARDIAN® INVISIBLE RETAINERS

GUARDIAN® Retainers are processed from your scans, models, or impressions. From your records, we will produce a printed model for laboratory fabrication at no charge. GUARDIAN® Retainers are available in quantities 2,3, and 4 per arch and are fabricated all at one time. All GUARDIAN® Retainers are processed on printed plastic models for maximum accuracy. These models are returned to give your practice the option to fabricate replacement retainers as needed. Zendura™ is used as a standard for all GUARDIAN® Retainers.

GUARDIAN® ENHANCEMENTS

Remove Brackets

We routinely remove all brackets on GUARDIAN® cases unless requested otherwise. Using the latest in advanced software, all fixed appliances are removed, free of charge, at the scanning stage to be reflected in the printed model. We recommend that all wires and bands are removed prior to submitting scans or impressions. If wires or bands are present, it will result in additional charges to remove them.

Add Fixed Lingual Retainer to a GUARDIAN® Series

Specialty Appliances will build a Fixed Lingual Retainer (FLR) design of your choice and process the GUARDIAN® Retainers to fit over the FLR. At the appointment when the braces are removed, the FLR is delivered with the GUARDIAN® Retainers fitting precisely over the FLR. Several designs of Fixed Lingual Retainers are available. Please use the diagram to indicate your wire and pad preferences.

Add Pontics to GUARDIAN® Retainers

Pontics can be added to any GUARDIAN® Retainer. We can also add pontics to any FLR design as a semi-permanent solution for patients too young for implants or restorations.

We encourage practices to contact our customer service representatives for additional information regarding these options.